

Bacterial Transformation

1.0
BETA Instructions and background notes

Transformation: a clue to life's mystery

Natural 'transformation' of bacteria was first described by the British microbiologist Fred Griffith in 1928. He showed that one strain of *Streptococcus* (then known as *Pneumococcus*) could be converted into another by an unknown, non-living material. The nature of Griffiths's 'transforming principle' remained a mystery for the next 15 years.

Oswald Avery and his colleagues in the USA carried out meticulous experiments for more than a decade to reveal the identity of the mystery material. Unfortunately, even though their findings (published in 1944) look conclusive to a modern observer, many scientists at the time argued against the Avery team's conclusion that DNA was the molecule of inheritance. Its molecular structure was thought to be too simple to carry the genetic message.

It was not until 1952 — just a year before the publication

of Watson and Crick's famous 'double helix' letter to *Nature* — that the majority of scientists were convinced that DNA was the primary genetic material. It was then that Alfred Hershey and Martha Chase showed, using radiolabelled DNA and protein, that the infective component of a bacteriophage was nucleic acid, and not the structurally more complex protein.

By the mid-1970s, transformation, which had provided a vital clue to the molecular nature of the gene, had become a key tool for the genetic modification of living things.

Genetic modification is central to many developments in modern biotechnology, and has been largely responsible for its evolution from a craft-based industry allied to brewing and baking to a major influence on human health and agricultural production in this millennium.

What you will learn

Educational aims

The following introductory practical procedure highlights several important ideas and provides:

- a practical demonstration that DNA is the genetic material;
- practical experience of one of the key techniques used in synthetic biology (bacterial 'transformation');
- an opportunity to learn, understand and carry out basic microbiological techniques;
- a concrete context for discussion of some of the ethical, social and safety issues associated with genetic modification;
- an opportunity to plan and carry out some more open-ended practical investigations (see page 20).

Summary of the practical task

You will transform a laboratory strain of *Escherichia coli* with plasmid DNA. The plasmid contains a gene encoding a green fluorescent protein (GFP) from the jellyfish *Aequorea victoria*. GFP glows brightly when illuminated with ultraviolet light, so it acts as a 'reporter' to confirm that the bacteria have indeed been transformed.

Bacterial transformation is a relatively inefficient process, and only a small proportion of the *E. coli* cells take up plasmid DNA. An antibiotic is therefore needed in the agar medium to prevent the growth of untransformed bacteria which would, if they were able to grow, outnumber the

transformed cells. The transformed cells can grow in the presence of the antibiotic, because a second gene on the introduced plasmid confers on its hosts resistance to the antibiotic.

Also on the plasmid is an origin of replication that allows the plasmid to be duplicated within the cells and a control region that 'switches on' the GFP gene. The switch is triggered by a chemical called IPTG, so in addition to the antibiotic, IPTG is incorporated into the agar medium.

The procedure used to transform the cells is known as the rapid colony method, and it uses a transformation buffer known as 'transformation and storage solution' (TSS). This method was first described by Chung *et al* in 1989 (see reference on page 8). It is not very efficient, but it is easy to carry out and does not require specialist equipment.

'Controls'

Normally, when undertaking a bacterial transformation such as this, one would carry out several 'control' treatments. These would include, for example, a 'transformation' without plasmid DNA. Several plates and types of agar media would be needed to perform all the necessary tests.

You may like to consider what sort of control treatments would be required, and carry out one or more of them if time and resources permit. In addition, you may like to investigate the effect of changing several aspects of the practical procedure: some ideas are given on page 20 of this guide.

Equipment and materials

Each person or working group will need

Straight from the kit

- a copy of these instructions
- a sterile single use spreader
- a sterile single use 5 μ L inoculation loop
- access to a UV LED torch (after the plates have been incubated)

Prepared in advance

(see pages 4, 5 and 6 of this guide for instructions)

- access to a stock culture of *E. coli*, K-12 strain TG2, **prepared no more than 48 hours in advance** (to be shared by group)
- 200 μ L of plasmid DNA in transformation buffer (TSS) dispensed into a microcentrifuge tube, on ice*.
- a Petri dish containing LB agar, kanamycin and IPTG*.


Not in the kit: supplied by you

- a small insulated cup of crushed ice
- a micropipette (e.g., 40–200 μ L)
- sterile pipette tips (for the micropipette)
- discard jar of freshly-diluted 1% (w/v) *Virkon*® disinfectant (for disposal of contaminated waste)*
- a permanent marker pen
- access to an incubator set at 37 °C
- a Bunsen burner

* Safety data sheets are provided for these items.

IMPORTANT

You must wear a lab coat and follow good microbiology laboratory practice while carrying out the practical work (see pages 13–15 of this guide).



You will also need:

- a copy of these instructions
- access to a fresh stock culture plate of *E. coli* (strain TG2)
- access to an incubator set at 37 °C
- access to antibacterial soap and paper towels

Preparing the materials

Plasmid DNA and transformation buffer

The kit contains *concentrated plasmid DNA* and *transformation buffer* (TSS). The concentrated plasmid DNA must be diluted in the transformation buffer and dispensed into tubes for individual students or groups to use shortly before the practical session.

We have provided the plasmid DNA and transformation buffer in *two* sets of tubes. This is so that, should you so desire, you can use half of it for one group and half for another group at a later date.

The plasmid DNA is provided in 2 mL screw-capped tubes. Each tube contains 100 μ L of highly concentrated plasmid solution. Each tube of transformation buffer contains 2 mL of liquid.

The diluted plasmid DNA should be dispensed *no more than 24 hours* before the practical session. Once this has been done, the tubes of diluted plasmid must be stored at -18 to -20°C , in a freezer, until they are required. This is because once it has been mixed with the transformation buffer, the plasmid DNA may deteriorate if it is kept at room temperature for more than a few hours.

To dilute the plasmid DNA

1. Take one of the tubes of concentrated plasmid DNA and one of the tubes of transformation buffer from the freezer. Allow the liquids to thaw, on the bench at room temperature, for 10–15 minutes. (If you leave

the unmixed solutions for several hours on the bench, even in a warm room, they will come to no harm — it's only once they are mixed that the plasmid DNA may deteriorate.)


2. Tap the closed tube of concentrated plasmid solution firmly on the bench a few times to return all of the liquid to the bottom of the tube.
3. Use a micropipette with a sterile tip to draw up *all* of the plasmid concentrate and add it to the tube of transformation buffer.
4. Cap the tube of liquid tightly and invert it several times to mix the contents.

To dispense the diluted plasmid ready for use

Each student or working group will need $\sim 200\mu\text{L}$ of the diluted plasmid DNA. For your convenience, we have provided a foam block in the kit for holding the tubes.

1. Use a micropipette with a sterile tip to dispense 200 μL of the diluted plasmid solution into each of the microcentrifuge tubes.
2. Close the tubes firmly after dispensing the liquid.
3. Tap each tube firmly on the bench to ensure that all the liquid lies at the bottom of the tube, then store the tubes in a freezer at -18 to -20°C until they are required.

IMPORTANT: If you store the frozen, diluted plasmid solution for more than 24 hours, it may deteriorate.



Growth media

This kit contains two sorts of agar growth media, in sachets: LB agar without added antibiotic and LB agar with added kanamycin and IPTG. Unlike conventional media, these must be prepared in a *microwave oven*. They *must not* be prepared by autoclaving as this would destroy the kanamycin in the medium. **IMPORTANT: Used plates must be disposed of in the normal way by autoclaving them.**

Each sachet contains sufficient material to make 200 mL of agar medium. The instructions below should be followed carefully to ensure good results.

Preparation of the sterile growth media

You will need:

- a microwave oven
 - heat-proof gloves
 - 2 × 500 mL borosilicate glass flasks, beakers or bottles
 - 1 × 1 L borosilicate glass flask or bottle
 - plastic film to cover the glass containers
 - 400 mL of distilled or deionised water
 - a water bath at 50 °C
1. Take the sachet labelled 'LB agar base'. Empty the contents into a clean 500 mL borosilicate glass bottle, flask or beaker.
 2. Add 200 mL of deionised or distilled water. *If you are using a beaker, cover it with plastic film, punctured once or twice.*
 3. Heat the liquid in a microwave oven on a MEDIUM power setting until bubbles start to appear (about 2–3 minutes). **IMPORTANT: You must watch the liquid constantly as it is heating to ensure that it does not boil over.**
 4. While wearing heat-proof gloves, take the container from the oven and swirl the liquid gently to mix. **CAUTION: Any solution heated in a microwave oven may become superheated and boil vigorously and suddenly when moved or touched. Take great care when handling the containers and ensure that you wear heat-proof gloves.**
 5. Reheat the liquid for 30 seconds at MEDIUM power and swirl gently as before. Do not overboil.
 6. Repeat step 5 if necessary until the powder is completely dissolved (the liquid should be clear, not turbid).
 7. Repeat steps 1–6 with the sachet of 'LB/X-Gal/Kanamycin/IPTG' agar.
 8. You should now have two containers of liquid media, which should be covered with clean plastic film or foil. Let the media cool to 50 °C (until the containers can be held comfortably in your hands).
 9. The media must be kept molten in a water bath at 50 °C until you are ready to pour the plates. **Remember to label the two containers in some way so that you can tell them apart.**
 10. While the agar medium is cooling, prepare the Petri dishes for pouring. **Do not pour the plates if the medium is hotter than 50 °C, as this will result in excessive condensation.**

Preparation of the agar plates

Each student or working group will require one Petri dish containing 12–15 mL of sterile kanamycin-containing agar.

In addition, you will need to prepare two or more plates of plain LB agar (without kanamycin) on which to grow cultures of untransformed bacteria that can be used as a source of bacteria by the entire class.

1. Open the bag of sterile Petri dishes. Cut the end of the plastic bag carefully so that it can be re-used to store the poured plates. Spread the Petri dishes out on the bench, unopened, ready to pour the agar.
2. Swirl the **LB agar** flask carefully to mix the agar. (When the mixture was prepared, the agar may have sunk to the bottom of the flask.)
3. When the agar has cooled to 50 °C, lift the lid of a Petri dish just enough to pour in some agar. Do not put the Petri dish lid down on the bench. Quickly add enough liquid to cover the bottom of the plate (you will need between 12 and 15 mL per Petri dish). Replace the lid immediately and tilt the plate to spread the agar.
4. Pour a second plate in the same way. Label both plates so that you know that they contain plain LB agar. *Note: If you wish to pour a few more stock plates of plain LB agar for class use, you may do so, but you may have to supply extra Petri dishes.*
5. Now mix the remaining LB agar with the **kanamycin-containing agar** in a larger, clean container. Swirl the container to mix the contents. *This is to dilute the kanamycin so that transformed bacteria will grow in its presence. The exact concentration of the kanamycin is not important.*
6. Continue pouring the rest of the plates, using this mixture. Flame the mouth of the container occasionally to maintain sterility.
7. If necessary, remove bubbles from the surface of the poured agar plates by briefly touching the surface with a Bunsen burner flame while the agar is still molten.
8. Allow the agar to solidify, undisturbed (this takes about 15 minutes).
9. Mark the plates so that you know they contain kanamycin. *Hint: an easy way to do this is to stack the plates, then to run a felt pen down the side of the stack, marking the edge of each Petri dish.*
10. Stack the plates, with the medium uppermost, in the original plastic sleeve for storage.
11. Surplus media may be used to produce extra plates or disposed of. **IMPORTANT: unused media MUST be autoclaved before disposal, to destroy the kanamycin it contains.**

Ideally, Petri dishes of agar should be poured *at least* 24 hours before the practical session, so that any contamination can be identified and to allow the agar to dry sufficiently. Plates that are contaminated **MUST NOT** be used. They should be disposed of by autoclaving.

Plates may be stored, inverted, (agar uppermost) for up to two weeks in a fridge at 3–5 °C.

IMPORTANT: Plates should be removed from the fridge a few hours before the practical session, to allow any condensation in the Petri dishes to evaporate.

Preparation and maintenance of the *E. coli* culture

The kit includes a slope culture of *Escherichia coli* TG2. This is the **ONLY** culture that should be used for this practical investigation. TG2 should be grown on LB agar as this contains all the nutrients the bacteria need.


Storage of the culture

Microbial slope cultures should generally be stored at *room temperature*, in the dark. They can be kept like this for up to 12 weeks, after which the bacteria should be subcultured onto a new slope of LB agar to ensure that you maintain a culture of viable cells. Slope cultures should **NOT** be stored in a fridge.

Streaking out

When you receive the culture, it is good practice to streak it onto a plate, progressively diluting the cells by spreading them out so that there are individual colonies on part of the plate. Stock plates should be prepared by taking cells from these colonies.

There are two common patterns for streaking plates: these are shown in the diagrams below. The (wire) loop should be flamed then allowed to cool between each set of streaks.



Prepare the plates one or two days in advance


Fresh stock cultures should be prepared from the slope by streaking the bacteria onto plates. Once you have prepared these plates they should ideally be used *within 48 hours*. To achieve good results, it is essential that a reasonable mass of cells is taken from the stock plate, taking care not to scrape any agar from the plate.

[If you keep the cultures for five days, you will still obtain results, but there will be far fewer transformed cells. If you keep the plates for a week, you *may* still obtain transformed cells, but typically there will be only three or four transformed colonies per plate.]

Subculturing for storage

Normally, you would subculture bacteria from a slope every eight to twelve weeks. In this way, you can usually maintain bacterial cultures for many months.

The strains used for genetic modification are often enfeebled, however, (e.g., lacking the enzymes needed for DNA repair if they are exposed to ultraviolet light). TG2 cells can also lose their naturally-occurring plasmid in storage. There is a risk that, after repeated sub-culturing, strains may accumulate undesirable mutations that will affect the success of the practical procedure. Unless you have the facilities for maintaining stock cultures at -80°C , it is better to obtain a fresh slope culture of TG2 when you wish to carry out this practical work.



Above: Two patterns which are often used when streaking out bacteria on plates.

Antibiotic resistance

The need for a selectable marker gene

The transformation process is very inefficient and only a small proportion of the bacteria treated will take up the novel plasmid. A means of selecting those cells that have been transformed is therefore needed.

Consequently, a gene is included on the plasmid which confers resistance to the antibiotic kanamycin, so that when kanamycin is included in the agar medium, only those bacterial cells that carry the plasmid will grow.

Resistance concerns

Resistance to the effects of many antibiotics is now widespread in several species of disease-causing microbes. It is important to appreciate how such resistance is transferred and selected for, and the special steps that have been taken in this particular kit to ensure that we do not contribute to this increasingly serious problem.

Resistance genes have evolved to give the bacteria that possess them the ability to thrive in environments containing antibiotics secreted by other microorganisms. Usually these genes encode enzymes that inactivate specific antibiotics or prevent them from working in some way. The widespread and often indiscriminate use of antibiotics in the treatment of disease has favoured resistant organisms, as the bacteria that are susceptible are simply killed off.

Natural transfer

Resistance genes are often carried on plasmids, which can pass from one bacterial cell to another of the same or a related species by a natural 'mating' process called conjugation. During conjugation, a tube or pilus is formed between adjacent cells, through which the plasmid passes.

The genes required for the formation of the pilus are also carried on a plasmid (an 'F' or fertility plasmid). The host strain provided in this kit (TG2) has an F plasmid, but this has been deliberately disabled so that no pilus is formed and consequently plasmids cannot transfer between cells.

Missing genes

For a plasmid to travel through a pilus, two additional requirements must be met. The plasmid must possess a gene encoding a mobility protein (*mob*) and have a *nic* site.

The mobility protein nicks the plasmid at the *nic* site, attaches to it there and conducts the plasmid through the pilus. The plasmid used in the transformation (*pNCBE-kan-GFP*) has neither a *nic* site nor the *mob* gene.

This ensures that once it has been introduced into a bacterial cell by artificial means (transformation) the novel plasmid cannot transfer into other bacterial cells.

As a further precaution, the plasmid DNA is non-methylated. That is, it is not protected from naturally-occurring restriction enzymes by methyl ($-CH_3$) groups. Therefore it would be degraded by restriction enzymes if it was to enter a wild-type bacterial cell.

The host strain TG2 is deliberately enfeebled, and has been selected for its non-pathogenic nature, its inability to survive outside the laboratory or to colonise the mammalian gut. It has several features that have been deliberately introduced into the genome of the bacterium to further assure biological safety (see *Bacterial strains and their genes*, pages 9–10 of this guide).

Physical containment

In addition to the numerous biological containment measures described above, the kit protocol requires that good microbiology laboratory practice is followed to ensure that the microorganisms are physically contained during the investigation and destroyed afterwards. Good microbiology laboratory practice is described on pages 13–15 of this document.

These methods of physical and biological containment have been adopted to make this educational protocol as safe as possible.

How kanamycin kills bacteria

Kanamycin kills bacteria by stopping protein synthesis in their cells. It does this by binding irreversibly to the 30S subunit of the bacterial ribosomes (eukaryote ribosomes have a different structure, so humans *etc.* are not affected). The kanamycin/ribosome complex initiates protein synthesis by binding to mRNA and the first tRNA. However, the second tRNA cannot bind, and the mRNA/ribosome complex dissociates.

Unlike several other antibiotics (*e.g.*, ampicillin) kanamycin kills all cells, rather than just those that are actively growing. Bacteria transformed with *pNCBE-kan-GFP* therefore require a short 'recovery period' before they are plated out onto kanamycin-containing plates. This allows the enzyme conferring resistance to be expressed.

How the kanamycin resistance works

Resistance to the effects of kanamycin is provided by the gene *kan^R* on the *pNCBE-kan-GFP* plasmid. This gene encodes an aminoglycoside 3'-phosphotransferase (APH).

APH catalyses the transfer of a phosphate group from ATP to a hydroxyl group of kanamycin. The phosphorylated antibiotic is unable to bind to the bacterial ribosome, so the antibiotic is inactivated. APH is relatively unstable, and is inactivated readily by increased temperatures or pH

changes. Its requirement for ATP means that this enzyme can only function in environments where that compound is abundant (e.g., inside cells) [For more information, see: www.rcsb.org/pdb/101/motm.do?momID=146]

Why use kanamycin?

While it would be easier to use a resistance marker that did not require a recovery period, there are several compelling reasons for using kanamycin and kanamycin resistance:

- unlike ampicillin and many other antibiotics, kanamycin is very seldom used to treat human disease, having been superseded by other drugs;
- it is needed in small amounts in culture plates (about a quarter of the concentration normally used for ampicillin);
- unlike ampicillin, kanamycin is not absorbed by the gut (in clinical use, it has to be injected). Therefore the safety hazard posed by accidental ingestion is reduced;
- kanamycin is relatively stable, so that plates containing it can be conveniently prepared well before a lesson;
- ampicillin resistance genes (β -lactamases) often confer resistance to other related antibiotics whereas the *kan^r*

gene affects a lesser range of antibiotics of limited therapeutic use;

- for several reasons, the use of kanamycin resistance markers is now widely accepted as safe, whereas scientists disagree about the wisdom of using ampicillin resistance markers.

In fact, the antibiotic resistance marker that is incorporated into *pNCBE-kan-GFP* might be better described as a kanamycin *tolerance* marker, as it does not confer full resistance to kanamycin to bacteria that possess it.

As stated above, the *kan^r* gene encodes an enzyme that modifies kanamycin, preventing it from working. This enzyme does not, however, alter all of the kanamycin in the medium, so high doses of kanamycin will prevent the growth even of the transformed bacteria. This is the reason that the LB/kanamycin/IPTG agar medium has to be 'diluted' with plain LB agar before use, to reduce the kanamycin concentration.

It is important, however, that all media that contain kanamycin are autoclaved before disposal to destroy the antibiotic, whether or not the media has been used to grow bacteria.

Further information

General reading

A glow in the dark by Vincent Pieribone and David F. Gruber (2005) The Belknap Press of Harvard University Press. ISBN: 978 0 674 02413 7. *Authoritative and beautifully-illustrated small book on the discovery and application of GFP.*

Glowing genes: A revolution in biotechnology by Mark Zimmer (2005) Prometheus Books. ISBN: 978 1591022534. *Popular but slightly error-prone account of the discovery and use of GFP.*

The transforming principle. Discovering that genes are made of DNA by Maclyn McCarty (1986) W. W. Norton and Company. ISBN: 978 0393304503. *Biographical account of the work of Avery, McCleod and McCarty by one of the participants.*

Web sites

In 2008, the Nobel Prize in Chemistry was awarded to Osamu Shimomura, Martin Chalfie and Roger Tsien for the discovery and development of GFP. More information can be found at the Nobel Prize web site: www.nobelprize.org/nobel_prizes/chemistry/laureates/2008/

Information about GFP and a template for a paper model can be obtained from the Protein Data Bank: www.rcsb.org/pdb/101/static101.do?p=education_discussion/educational_resources/GFP_activity.html

Bristol University produces a similar series which features an article on GFP written by Timothy King and Paul May: www.chm.bris.ac.uk/motm/GFP/GFP.htm

Molecular structure data

The computer-generated image on the cover of this booklet was created using structural data from the Protein Data Bank: www.rcsb.org/pdb

The image shows a computer model of a molecule of green fluorescent protein, using data from: Ormo, M., *et al* (1996) Crystal structure of the *Aequorea victoria* green fluorescent protein *Science* **273**, 1392–1395 [Protein Data Bank ID: 1EMA].

The software used to produce this image was UCSF Chimera, which can be obtained free-of-charge from: www.cgl.ucsf.edu/chimera/

Genetic modification method and microbiology safety

A guide to the genetically modified organisms (contained use) regulations 2000. Health and Safety Executive (2000) The Stationery Office, London. ISBN: 978 0717617586. *This official document, which is aimed principally at academic researchers, can be downloaded from the HSE's web site: www.hse.gov.uk/biosafety/gmo/*

Chung, C.T., Niemela, S. and Miller, R.H. (1989) One-step preparation of competent *Escherichia coli*: Transformation and storage of bacterial cells in the same solution. *Proc. Natl. Acad. Sci. USA*. **86**: 2172–2175. *This is the method for transforming cells used in this kit.*

Bacterial strains and their genes

Bacterial genotypes

The genotype of a bacterial strain is usually given by listing all of the genes that are known to differ from the wild type. This is done using a system proposed by Milislav Demerec and his colleagues in 1966. The main features of this system are:

- Each genetic locus in the wild type is designated by a three-letter, lower-case symbol, which is written in italics *e.g.*, *pro* is the symbol for the gene determining metabolism of the amino acid proline;
- Different loci are distinguished from one another by adding a capital letter after the symbol *e.g.*, *argA*; *argB*. If the exact locus at which the change (mutation) has occurred is not known, a hyphen is used instead of a capital letter *e.g.*, *ara-*;
- If it is known, the site at which the mutation has occurred is shown by a number after the locus letter *e.g.*, *hisA38*.
- Δ (Greek symbol delta) indicates that the genes following it have been *deleted*.
- Conversely, a double set of colons, ::, indicates that the genes following it have been *inserted*.

Phenotypic traits are described in words, or by abbreviations which are explained when they are first used. These abbreviations are clearly distinguished from those referring to the genotype *e.g.*, Arg- is a phenotype, showing that arginine is required, but *argA* is a specific locus at which a mutation has occurred.

Under this naming system, individual strains are designated by serial numbers, that are decided by the laboratories that have isolated those strains. These numbers are not italicised. For instance, *Escherichia coli* K-12 CSH50 is strain 50 from Cold Spring Harbor Laboratory, a well-known research establishment in the USA.

Demerec and his co-workers also devised a similar system for naming the genes on plasmids. To avoid confusion the plasmid genotype is listed within square brackets.

Escherichia coli TG2

For this practical investigation, we have supplied a strain of *E. coli* called TG2. This strain is suitable for transformation by plasmids and grows well on broth and agar plates which are supplemented with all of the amino acids the bacterium needs (LB agar contains all of these amino acids).

TG2 is a K-12 strain. Unlike the wild type, K-12 strains of *E. coli* are unable to inhabit the mammalian gut. This strain's origins can be traced back to work in the USA in 1922. Biochemical and genetic studies by Edward Tatum in the 1940s made the strain popular with researchers, and after many millions of generations of laboratory cultivation, it is

now known to have undergone significant changes. These have altered the lipopolysaccharides that compose the outer membrane of the bacterial cell, so that it can no longer infect mammals (the cells lack the 'O' antigen which is required for infection). This inability to colonise the gut makes them particularly safe for laboratory use.

TG2 is one of many mutant strains of *E. coli* that has been developed specially for use in genetic modification. It carries the plasmid of its own and can be transformed efficiently by certain plasmids or bacteriophages (for instance, several special features have been introduced into the bacterium to assist its transformation with the bacteriophage M13).

Like most cloning strains, compared to the wild-type *E. coli*, TG2 is severely weakened and it would find it difficult to thrive outside the laboratory. For example, although it can grow on glucose, unlike the wild type *E. coli*, it is unable to use lactose as an energy source.

Genotype of TG2

The genes of this strain that are of interest to molecular biologists are usually represented by the following gene symbols:

Bacterial chromosome

supE, *hsd* Δ 5, *thi*, Δ (*srl-recA*)306::Th10(*tet*^r), Δ (*lac-proAB*)

Bacterial (F') plasmid

F' [*traD*36, *proAB*⁺, *lacI*^q, *lacZ* Δ M15]

A. Bacterial chromosome

supE (also known as *glnV*)

This mutation is a safety feature left over from the early days of genetic modification.

The triplet UAG is usually recognised as a 'STOP' codon. Bacteria with the *supE* mutation have mutant tRNA, however, and instead of stopping at UAG during protein synthesis they instead insert glutamine into the amino acid chain.

Genetic modification techniques sometimes involve the use of M13 bacteriophage (a virus that infects bacteria) to transfer genes into cells. In the past, these viruses were required by safety authorities to have UAG ('STOP') codons inserted into several important genes within the virus genome. When introduced into bacterial hosts with the *supE* mutation, such as TG2, the 'STOP' codons in the viral DNA were ignored, allowing the novel genes to be expressed. If, however, the viral DNA accidentally transferred into wild-type bacteria (that is, it 'escaped' into the environment) the numerous stop codons would prevent the viral genes from being translated into proteins.

hsdΔ5

Normally, *E. coli* produces restriction enzymes to cut up any foreign DNA that enters the cell. The bacterium's own DNA is protected from these restriction enzymes because it is methylated (methyl, $-CH_3$, groups are added to the DNA). The *hsd* genes encode enzymes responsible for both restriction and methylation of DNA.

hsdΔ5 is a mutation in the system of methylation and restriction that allows non-methylated plasmid DNA to be introduced into the bacterium without risk of it being degraded. The use of non-methylated plasmid DNA for genetic modification is another safety precaution, because should the plasmid enter a wild-type bacterium, without the protection that methylation affords, it will be recognised as foreign and will be broken down.

thi

This indicates that TG2 requires thiamine (vitamin B₁) for growth. Vitamin B₁ is present in LB and nutrient agar.

$\Delta(srl-recA)306$

$\Delta(srl-recA)306$ is a deletion from the bacterial genome stretching from the *srl* genes through to the *recA* gene. The *srl* genes control the metabolism of sorbitol, and consequently TG2 is unable to use sorbitol as an energy source.

Numerous *rec* (recombination) genes are found in *E. coli*. The proteins they encode control recombination of the bacterium's DNA. If present, recombination proteins can also rearrange the DNA in a plasmid or bacteriophage (such as M13 or lambda) used to genetically modify the bacterium. Deletion of the *recA* gene prevents recombination of introduced DNA.

The number '306' after the brackets is simply to distinguish this deletion from the 305 other deletions that geneticists had previously made to the bacterium's genome.

::Tn10(*tet^r*)

A gene encoding resistance to the antibiotic tetracycline, *tet^r*, has been inserted. Hence TG2 is able to grow on media containing the antibiotic tetracycline. Tn10 indicates that the *tet^r* gene comes from the Tn10 transposon.

$\Delta(lac-proAB)$

The genes from *lac* to *proAB* have been deleted from the bacterial genome.

The *lac* gene encodes the lac operon, which is needed to trigger the production of β -galactosidase (lactase), which breaks down lactose. Consequently, unlike wild-type *E. coli*, TG2 is unable to metabolise lactose.

proAB symbolises two genes that bacteria need to synthesise the amino acid proline. In the absence of these genes, bacteria would normally require proline in their growth medium. Note, however, that TG2's F' plasmid includes the *proAB* genes, so as long as TG2 retains its own plasmid it will have the ability to make proline. [This transfer of the *proAB* genes from the bacterial chromosome to its plasmid has been done deliberately to make it easy to maintain cells with the bacterial plasmid in culture — see below.]

B. Bacterial (F') plasmid

TG2 has an F (fertility) plasmid. An F plasmid usually carries, amongst other genes, those encoding the proteins responsible for making the 'sex pilli' that permit genetic material to pass between cells.

The prime (') after the F indicates that the plasmid has also picked up some genes from the bacterium's chromosome.

traD36

This is a mutation in one of the genes involved in bacterial conjugation. The mutation has been introduced as a safety feature so that TG2 cannot naturally pass genetic material to other bacteria by conjugation. Therefore if TG2 is genetically modified with plasmid DNA, this novel genetic material cannot be passed to other bacteria.

The remaining genes: *proAB*, *lacI^q* and *lacZΔM15*, have been deleted from the bacterial chromosome and transferred to the F' plasmid. This has been done to ensure that the F' plasmid remains in the culture, to exercise greater control of the expression of introduced genes in genetically-modified bacteria and to help identify any bacterial cells that have been genetically modified.

***proAB*⁺**

These two genes return the ability to synthesise proline to TG2. F plasmids are sometimes lost from bacteria after repeated subculturing. By growing stock cultures of the bacteria on minimal medium (that contains no proline), scientists can ensure that only those bacteria that retain the F' plasmid will survive.

lacI^q

The *lacI* gene encodes the repressor of the *lac* operon (the *lac* operon is the genetic 'switch' that activates the production of β -galactosidase).


lacI^q is a mutant form of the gene that produces about ten times more repressor than the wild type, and it does this continuously (the 'q' indicates that the mutation is constitutive: a constitutive gene is one that is transcribed continually). This ensures that any genes introduced into the bacteria that are under control of the *lac* operon are *not* expressed unless the inducer, IPTG, is present in the growth medium. This mutation can therefore be thought of as a strong 'off' switch.

lacZΔM15

lacZ is the gene encoding β -galactosidase (lactase), which would normally enable the bacterium to metabolise lactose. The mutant gene on the plasmid (*lacZΔM15*) lacks codons 11 to 41. These codons specify part of the protein called the α -peptide, the absence of which leads to the production of a non-functional β -galactosidase.

To make β -galactosidase, therefore, the missing part of the gene (encoded by *lacZ'* and sometimes called α -complement) has to be provided. This is often done by including the necessary fragment of the gene on plasmids or forms of the M13 bacteriophage used for genetic modification.

Outline plasmid map



Origin of replication

The origin of replication allows the plasmid to be reproduced within the bacterium, producing between 50 and 700 copies of the plasmid within the cell (depending upon the conditions of growth).

Control region

This is a 'genetic switch' which activates transcription of the GFP gene. IPTG in the growth medium mimics the sugar lactose, activating the 'switch'.

Green fluorescent protein (avGFP)

This marker gene is derived from the jellyfish *Aequorea victoria*. The green fluorescent protein (avGFP) it encodes glows brightly in near-ultraviolet light of ~385 nm.

Kanamycin resistance (*Kan^R*)

The gene *Kan^R* encodes aminoglycoside 3'-phosphotransferase (APH). This enzyme alters the antibiotic kanamycin, preventing it from working.

Safety and genetic modification



Contained Use

All practical work that involves the production or use of genetically-modified organisms (GMOs) is strictly regulated by law throughout the European Union. There are two principal sets of EU regulations (Directives) governing genetic modification. Laws in the United Kingdom and elsewhere within the EU are enacted to comply with these Directives. One Directive [Directive 2009/41/EC] covers 'Contained Use' *e.g.*, work in a laboratory; the other [Directive 2001/18/EC] covers the 'Deliberate Release' of GMOs into the environment *e.g.*, field trials of genetically-modified crops.

In general, anyone carrying out work with GMOs must do so only on premises that have been registered with the relevant authority. In the UK, this is principally the Health and Safety Executive (HSE). The organisation, such as a university or research facility, under whose auspices the work is to be done must also set up a local expert safety committee and procedures to oversee and control the work with GMOs. These stringent requirements mean that in the European Union almost all work with GMOs is restricted to universities and research institutions. There is, however, a limited amount of practical work that is exempt from the Contained Use Regulations.

'Self-cloning'

The practical procedure described in this kit is known technically as 'self-cloning'. Here, 'cloning' means making copies of DNA within an organism. Originally, the definition of self-cloning was restricted to taking DNA from one species and making copies of it (cloning it) in the *same* species — hence the term *self-cloning*. Later, this definition was widened slightly to include 'reporter' or 'marker' genes and control sequences, which might come from other species, provided these elements had an extended history of safe use.

Self-cloning using non-pathogenic microorganisms, such as the weakened laboratory strain *E. coli* provided with this kit, is *exempt* from the Contained Use regulations. The bacteria produced *are* covered by the Deliberate Release regulations, however, and it is therefore essential to ensure that an accidental 'release' of the organism into the environment does not occur. This is achieved in two ways: by physical and by biological containment.

Physical containment

The genetically-modified microorganisms (GMMOs) must be *physically contained* by good microbiology laboratory practice, including the destruction of the cultures after use. The relevant basic microbiology laboratory techniques are described on pages 13–15 of this document. These should be read carefully and followed by those undertaking the work.

Biological containment

The GMMOs are also *biologically contained*, by the selection of a suitable host strain and the careful construction of the plasmid DNA. So, for example, in the current practical procedure, the strain of *E. coli* lacks the ability to pass on the introduced DNA by the natural bacterial 'mating' process of conjugation, and the plasmid DNA is non-methylated so that if it *did* enter a wild-type bacterium, it would be degraded by that organism's own restriction enzymes. Further details of the genotype of the host strain and the plasmid construction are provided on pages 9–11 of this document.

Alterations to the procedure

It follows from what has been stated above that *no attempt* should be made to alter or add to the procedure described in this kit in a way that might bring those following it outside the umbrella of self-cloning and into the realm of Contained Use. If this was to be done without notifying the relevant authorities and following the other procedures that such work legally requires, users could place themselves and others at risk, and could ultimately be subject to legal action.

Further information

Additional information and guidance regarding health and safety relating to work with GMOs can be found on the HSE's web site: www.hse.gov.uk/biosafety/gmo/ and in the following publication: *A guide to the genetically modified organisms (contained use) regulations 2000*. Health and Safety Executive (2000) The Stationery Office, London. ISBN: 978 0717617586. This document can be downloaded from the HSE's web site.

Note that this publication is aimed primarily at research institutions and is concerned with technical procedures relating to 'Contained Use', often on a large scale or involving hazardous microorganisms.

Video demonstrations of basic microbiology laboratory techniques and other useful information can be found on the Society for General Microbiology's YouTube web site: www.youtube.com/user/SocGenMicrobiology/videos

Good microbiology laboratory practice

General precautions

- Any exposed cuts and abrasions should be protected with waterproof dressings before the practical work starts.
- There is *no need* to wear disposable gloves, except if a person has skin condition such as eczma or abrasions or cuts to the skin that cannot be covered with waterproof dressings (either because they are too large or awkward to cover or the person concerned is allergic to plasters).
- Everyone involved — lecturers, technicians and students should wash their hands *before and after* practical work.
- Laboratory doors and windows should be closed while practical work is in progress. This will reduce air movements and consequently the risk of accidental contamination of plates, *etc.*
- High standards of cleanliness must be maintained. Non-porous work surfaces should be used and they must be swabbed with an appropriate laboratory disinfectant before and after each practical session. (*Virkon*[®] is the disinfectant of choice for microbiology).
- No hand-to-mouth operations should occur (*e.g.*, chewing pencils, licking labels, mouth pipetting). Eating, drinking and smoking must not be allowed in the laboratory.
- Those carrying out the work should wear laboratory coats and, where necessary, eye protection.

Spills and breakages

Accidents involving cultures should be dealt with as follows:

- Disposable gloves should be worn.
- The broken container and/or spilt culture should be covered with paper towels soaked in disinfectant.
- After not less than 10 minutes, it must be cleared away using paper towels and a dustpan.
- The contaminated material must be submerged in a suitable disinfectant for 24 hours or placed in a microbiological disposal bag.
- The contaminated material must then be autoclaved before disposal. The dustpan should also be autoclaved or placed in a bucket of suitable disinfectant solution (*e.g.*, *Virkon*[®]) for 24 hours.
- Contaminated paper towels should be autoclaved.

Contamination of skin or clothing

As soon as possible, anyone affected should wash with antibacterial soap. Severely contaminated clothing should be placed in disinfectant before it is laundered.



Sources of microbes

All micro-organisms should be regarded as potentially harmful. However, the strain of *E. coli* used in this kit presents minimum risk given good microbiology laboratory practice. Other species of bacteria *must not* be used for this work, as this might contravene the regulations governing genetic modification (see page 12).

In general, stock slope cultures of bacteria should be kept in the dark at room temperature, not in a fridge. Slope cultures should be subcultured onto fresh nutrient agar every six weeks or so. You should not attempt to maintain the culture for an extended period, however, as mutations can occur in the storage conditions that are found in schools, and these may lead to the failure of the practical work. Cultures may also become contaminated with repeated sub-culturing. If in doubt, obtain a fresh culture.

Aseptic techniques

The aims of aseptic techniques are:

- To obtain and maintain pure cultures of microorganisms;
- To make working with microorganisms safer.

A 'pure culture' contains only one species of microorganism, whereas a 'mixed culture' contains two or more species.

Contamination of cultures is always a threat because microbes are found everywhere; on the skin, in the air, and on inanimate objects. To obtain a pure culture, sterile growth media and equipment must be used and contaminants must be excluded. These are the main principles of aseptic techniques.

It is unrealistic to expect inexperienced school students to be fully accomplished at aseptic techniques. Sterile, disposable items are therefore provided in this kit, so that the necessary procedures can be carried out as easily and safely as possible.

Growth media must be prepared as described on page 5 of this booklet. Sterile Petri dishes should be used. Lids must be kept on containers to prevent contamination.

Practical work should be carried out near a Bunsen burner flame. Rising air currents from the flame will carry away any microbes that could contaminate growth media and pure cultures.

When cultures are transferred, tops and lids of containers should not be removed for longer than necessary. After a lid has been taken from a bottle, it should be kept in the hand until it is put back on the bottle. This prevents contamination of the bench and the culture. After removal of the top, the neck of the culture bottle should be flamed briefly for 1–2 seconds. This will kill any microbes present there and produce convection currents which will help to prevent accidental contamination of the culture.

With practice, it is possible to hold a bottle containing the microbes in one hand and the loop or pipette in the other in such a way that the little finger is free to grip the bottle top against the lower part of the hand. (In this case, it is important that the bottle top should be loosened slightly before the inoculation loop is picked up.)

Obviously, unlike glassware, the sterile plastic spreaders, loops, microcentrifuge tubes and pipettes that are provided in this kit must *not* be flamed.

If, however, you use wire loops (*e.g.*, as replacements for the disposable loops in this kit) they should be heated until they are red hot along the entire length of the wire part. This should be done both before and after transfer of cultures takes place. Loops should be introduced slowly into the Bunsen burner flame to reduce sputtering and aerosol formation.

When the Bunsen is not in use it should be kept on the yellow flame, so that it can be seen. A blue flame about 5 cm high should be used for sterilising loops and flaming the necks of bottles.

Avoid contaminating the work area. Any non-disposable instruments should be sterilised immediately after use and used plastic pipettes, tubes and other plastic items should be placed directly into a jar of fresh disinfectant (*Virkon*®) solution, so that they are completely immersed.

Incubation of cultures

Label the Petri dish around the edge of the base. A name, date and the name of the organism used will allow the plate and its contents to be identified.

Bacterial cultures in Petri dishes should be incubated with the base uppermost, so that any condensation that forms falls into the lid and not onto the colonies.

Disposal

It is very important to dispose of all materials properly after use, especially any item that has been in contact with cultures of bacteria. All non-disposable containers used for storing and growing cultures must be autoclaved, then washed and rinsed as necessary, before re-use.

There should be a discard jar of fresh disinfectant (*Virkon*®) near each work area. Disposable plastic pipettes, loops, spreaders, tubes and any liquid from cultures should be put into the disinfectant pot immediately after use. After soaking for 24 hours, these materials should be autoclaved then disposed of in the normal waste.

Contaminated paper towels, cloths and plastic Petri dishes should be put into an autoclave bag and sterilised by autoclaving before placing in the normal waste.

Glassware that is not contaminated (*e.g.*, flasks used for making up media) can be washed normally. Broken glassware should be put in a waste bin reserved exclusively for that purpose. If any glassware is contaminated it must be autoclaved before disposal. Uncontaminated broken glassware can be thrown away immediately.

Autoclaving

Sterilisation is the complete destruction of all micro-organisms, including their spores.

All equipment should be sterilised before starting practical work so that there are no contaminants. Cultures and contaminated material should also be sterilised after use for safe disposal.

Autoclaving is the preferred method of sterilisation for culture media, aqueous solutions and discarded cultures. To comply with the regulations governing genetic modification, it is *essential* that any genetically-modified microorganisms are killed before disposal. Autoclaving is the most reliable method of achieving this.


The process uses high pressure steam, usually at 121 °C. Microbes are more readily killed by moist heat than dry heat as the steam denatures their proteins. Domestic pressure cookers can be used instead of autoclaves, but their small capacity can be a disadvantage when dealing with larger amounts of material.

Principles of autoclaving

Two factors are critical to the effectiveness of autoclaving. Firstly, all air must be removed from the autoclave. This ensures that high temperature steam comes into contact with the surfaces to be sterilised: if air is present the temperature at the same steam pressure will be lower.

The materials to be sterilised should be packed loosely so that the air can be driven off. Screw-capped bottles and jars should have their lids loosened slightly to allow air to escape and to prevent a dangerous build-up of pressure inside them.


Secondly, enough time must be given for heat to penetrate (by conduction) to the centre of media in Petri



dishes or other containers. The times for which media or apparatus must be held at various temperatures for sterilisation are shown in the graph and table below. Notice that just a small difference in temperature can result in a significant difference in the time required for sterilisation. It is also important that these temperatures are reached by *all* of the materials to be sterilised for the specified time e.g., the broth in the very centre of a flask.

Three factors contribute towards the duration of the autoclaving process:

- **penetration time** — the time taken for the centre-most part of the autoclave's contents to reach the required temperature;
- **holding time** — the minimum time in which, at a given temperature, all living organisms will be killed;
- **safety time** — half the holding time, included as a safety margin.



Domestic pressure cookers operate at 121°C. Thus the total sterilisation time might typically be: a penetration time of 5 minutes plus 15 minutes holding time (including a safety margin) — giving a total time of 20 minutes.

The effectiveness of an autoclave can be checked by using autoclave test strips (available from laboratory suppliers) which change colour if the process has worked properly (autoclave tape does not show this).

Use and care of autoclaves

The manufacturer's instructions should always be followed when using a pressure cooker or an autoclave. Particular care should be taken to ensure that there is enough water in the autoclave so that it does not boil dry during operation. Most domestic pressure cookers require at least 250 mL of water — larger autoclaves may need far greater volumes.

The use of distilled or deionised water in the autoclave will prevent the build-up of limescale. Autoclaves should be dried carefully before storage to prevent corrosion of the pressure vessel.

When the autoclave is used, before the exit valve is tightened, steam should be allowed to flow freely from the autoclave for about one minute to drive off all the air inside. After the autoclave cycle is complete, sufficient time must be allowed for the contents to cool and return to normal atmospheric pressure. The vessel or valves should not be opened whilst under pressure as this may cause scalding.

Premature release of the lid and the subsequent reduction in pressure will cause any liquid inside the autoclave to boil and spill from vessels.

Chemical sterilisation

Historically, many different chemicals have been used in microbiology to sterilise used equipment and work surfaces. Hypochlorites (such as bleaches) and clear phenolics should no longer be used, however. *Virkon*® can be safely used for most microbiology laboratory purposes. The solution remains pink while active, but if it turns colourless or pale, it is no longer effective and should be disposed of.

The manufacturer's and supplier's instructions should always be followed with care. A safety data sheet for *Virkon*® is provided with this kit. **IMPORTANT: Eye protection must be worn when handling concentrated solutions of disinfectant.**

Handling disposable sterile items


The plastic loops and spreaders provided in this kit are sterile. The packets should be opened with care to ensure that the tips of pipettes, loops and so on are not touched.

These items should all be placed in a discard container of disinfectant after use, where they should be left for 24 hours before autoclaving and disposal in the normal waste.

Sterilising microcentrifuge tubes

If necessary, microcentrifuge tubes can be sterilised by autoclaving before use — they will not melt. Place them in a beaker, cover it with aluminium foil held on with a strip of autoclave tape and autoclave at 121 °C for 15 minutes.

After use, microcentrifuge tubes should be placed, open, in *Virkon*® disinfectant solution, left for 24 hours, then autoclaved before disposal in the normal waste.



Rapid colony transformation

Using Transformation and storage solution (TSS)

Before you start, put on a lab coat and wash your hands with antibacterial soap.

1. Wipe down the work surface using disinfectant (*Virkon*®) solution and paper towels. Set out the work area with everything you need (see page 3), including a container of disinfectant solution for contaminated waste and a Bunsen burner with a yellow flame.
2. Place the plasmid DNA (in TSS buffer) on ice.
3. Use a sterile loop to scrape a visible mass of bacteria from the stock plate (but don't dig into the agar).
4. Put the loop of bacteria into the ice-cold plasmid DNA solution and rotate it vigorously to dislodge the bacteria and suspend them in the solution. Hold the tube almost horizontally as you do so. Once you have dislodged the bacteria from the loop, place the contaminated loop in the disinfectant solution.
5. Using a sterile tip on the micropipette, gently draw the liquid in the tube up and down a few times to ensure that the bacteria and the plasmid are well mixed. Dispose of the used tip into the disinfectant.
6. Close the tube tightly and put it on ice for 15 minutes.
7. During this time, the plasmid DNA will be taken up by the bacteria.
8. Put a sterile tip on a micropipette and use it to drop 150 μL of bacterial suspension onto the centre of a Petri dish of LB/Kanamycin/IPTG agar. Place the contaminated pipette tip and open tube of bacteria in the disinfectant.
9. Use a sterile spreader to coat the surface of the agar with the bacterial suspension. Rotate the plate as you do this, so that the culture is spread evenly over the plate. Once you have done this, place the contaminated spreader in the disinfectant.
10. Write your name, the date and the bacteria used on the edge of the base of the Petri dish.
11. Place the Petri dish, inverted (agar uppermost) in an incubator at 37 °C. Clean the work area with disinfectant and wash your hands with antibacterial soap.
12. Incubate the plates overnight — ideally for 24 hours — then examine the bacteria under ultraviolet light. Colonies expressing GFP should fluoresce brightly.






After the plates have been examined, dispose of all contaminated materials and cultures by autoclaving them.




Transformation efficiency

How good were you at transforming bacteria?

Transformation efficiency is expressed as the number of transformed colonies produced per microgram (μg) of plasmid DNA.

1. Count the number of transformed colonies on your plate.
Number of colonies = _____
2. Calculate the mass, in micrograms (μg), of plasmid DNA used.
(There were 0.5 micrograms (μg) of plasmid DNA in the 200 μL of solution you were given.)
Mass of plasmid DNA used = _____
3. Note the volume of bacterial cell suspension you spread on your plate.
Volume of suspension spread on plate = _____
4. Calculate the fraction of the cell suspension you spread on your plate (volume of suspension spread [3] \div total volume of suspension (200 μL)).
Fraction of suspension spread on plate = _____
5. Calculate the mass of plasmid DNA (in μg) contained in the cell suspension you spread on the plate (total mass of plasmid [2] \times fraction spread [4]).
Mass of plasmid DNA spread on plate = _____
6. Calculate the transformation efficiency (number of colonies [1] \div mass of plasmid spread [5]).
Transformation efficiency = _____





Thinking about ‘controls’

1. What would you expect to see on the following plates and why?

- a. Bacteria to which plasmid DNA had *not* been added, that were plated on LB agar containing antibiotic and IPTG.

- b. Bacteria to which plasmid DNA had been added, that were plated on LB agar containing antibiotic but no IPTG.

- c. Bacteria to which plasmid DNA had been added, that were plated on LB agar containing IPTG but no antibiotic.

2. Describe the experiments that would you need to carry out to test whether the antibiotic kanamycin makes transformed bacteria glow in UV light.

Further investigations

Although for safety reasons, there is a limit to the scope of extensions to this practical work that can be attempted, there are several variations on the technique described in this booklet which can be carried out. For example, you could investigate the effect of changing one or more of the following conditions:

- the age of the host cells used;
- the mass of plasmid DNA used;
- the concentration of kanamycin in the medium (note that the 'undiluted' kanamycin agar contains 50 mg of kanamycin per mL of medium);
- the duration of the recovery period;
- the temperature of the recovery period.

Because of the nature of the kanamycin resistance (see pages 7–8) changing the duration of the recovery period

will change the transformation efficiency dramatically. A recovery period of two hours, for example, will produce many more bacterial colonies than a 15-minute period.

Transformation efficiency is expressed as the number of transformed bacterial colonies produced per microgram (μg) of plasmid DNA. The worksheet on page 16 guides you through the necessary calculation.

Some of these investigations could be attempted as whole class exercises. For example, if different students or groups in the class are each provided with a different concentration of plasmid DNA, graphs of DNA mass vs the number of colonies, and DNA mass vs the transformation efficiency can be plotted.

Such graphs can be used to determine the point at which the transformation reaction is saturated (the point at which increasing the concentration of plasmid no longer results in greater transformation efficiency).

Plasmids and DNA sequence data

The sequences of the plasmids and other inserts needed for these practical tasks are available from *Figshare*:
www.figshare.com [Search for 'Practical synthetic biology']

The plasmids and inserts will be available from the National Centre for Biotechnology Education, University of Reading, UK (**www.ncbe.reading.ac.uk**).

For those unable to obtain the plasmids from the NCBE, they have also been deposited in *AddGene*:
<https://www.addgene.org> [Search for 'Bryk']

Updated versions of the written and other educational resources may be found on the UNIGEMS project website:
practicalsyntheticbiology.net/resources

Acknowledgements



The *pNCBEkanGFP* plasmid and this practical protocol were developed by Jarek Bryk at the National Centre for Biotechnology Education (NCBE), University of Reading, UK. This project (UNIGEMS) has received funding from the European Union's Seventh Framework Programme for research, technological development and demonstration under grant agreement No. 300038.